

Sonaer Part Number	Syringe Size ml	Syringe Material	Max Pressure
SON00S20	0.05	Glass With PTFE Plunger	400
SON00S21	0.10	Glass With PTFE Plunger	400
SON00S22	0.25	Glass With PTFE Plunger	400
SON00S23	0.50	Glass With PTFE Plunger	400
SON00S02	1.00	Polypropylene	N/R
SON00S11	1.00	Standard Glass	43.5
SON00S34	1.00	Stainless Steel	750
SON00S24	1.00	Glass With PTFE Plunger	400
SON00S03	3.00	Polypropylene	N/R
SON00S13	3.00	Standard Glass	43.5
SON00S38	3.00	Stainless Steel	750
SON00S25	2.50	Glass With PTFE Plunger	400
SON00S04	5.00	Polypropylene	N/R
SON00S14	5.00	Standard Glass	43.5
SON00S26	5.00	Glass With PTFE Plunger	400
SON00S42	5.00	Stainless Steel	750
SON00S05	10.00	Polypropylene	N/R
SON00S15	10.00	Standard Glass	29
SON00S27	10.00	Glass With PTFE Plunger	200
SON00S46	8.00	Stainless Steel	1500
SON00S06	20.00	Polypropylene	N/R
SON00S16	20.00	Standard Glass	29
SON00S28	25.00	Glass With PTFE Plunger	100
SON00S51	20.00	Stainless Steel	750
SON00S07	30.00	Polypropylene	N/R
SON00S17	30.00	Standard Glass	29
SON00S08	50.00	Polypropylene	N/R
SON00S18	50.00	Standard Glass	21.8
SON00S29	50.00	Glass With PTFE Plunger	100
SON00S56	50.00	Stainless Steel	750
SON00S09	60.00	Polypropylene	N/R
SON00S19	100.00	Standard Glass	21.8
SON00S30	100.00	Glass With PTFE Plunger	100
SON00S61	100.00	Stainless Steel	750
SON00S10	140.00	Polypropylene	N/R

<u>Syringe Pump Kit Description</u>	<u>Qty.</u>
25' FEP 1/8" OD x 1/16" ID TUBING	1
ETFE LEUR TO 1/4"-28 FEMALE ADAPTER	1
ETFE NUT FOR 1/8" TUBE	1
1/8" FERRULE ETFE	1
1/8 NPT TO 1/8" OD TUBE COMPRESSION FITTING	2
1/8" OD TUBING COMPRESSION FITTING	1
1/8" OD TUBING BLACK NYLON/NICKEL-PLATED BRASS PUSH TO CONNECT COUPLING	1
1/8" OD TUBING BLACK NYLON/NICKEL-PLATED BRASS PUSH TO CONNECT ELBOW	1

<u>Gear Pump Kit Description</u>	<u>Qty.</u>
25' FEP 1/8" OD x 1/16" ID TUBING	1
1/4"-28 ETFE UNION	1
ETFE NUT FOR 1/8" TUBE	1
1/8" FERRULE ETFE	1
1/8 NPT TO 1/8" OD TUBE COMPRESSION FITTING	2
1/8" OD TUBING COMPRESSION FITTING	1
1/8" OD TUBING BLACK NYLON/NICKEL-PLATED BRASS PUSH TO CONNECT COUPLING	1
1/8" OD TUBING BLACK NYLON/NICKEL-PLATED BRASS PUSH TO CONNECT ELBOW	1